

EphA10 Polyclonal Antibody

Cat No: HR1AP6712

For research use only

Overview

Product Name	EphA10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	110kDa
GeneID?Human?	EPHA10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EPHA10; Ephrin type-A receptor 10
Background	EPH receptor A10(EPHA10) Homo sapiens Ephrin receptors, the largest subfamily of receptor tyrosine kinases (RTKs), and their ephrin ligands are important mediators of cell-cell communication regulating cell attachment, shape, and mobility in neuronal and epithelial cells (Aasheim et al., 2005 [PubMed 15777695]). See MIM 179610 for additional background on Eph receptors and ephrins.[supplied by OMIM, Mar 2008],