

Relaxin Receptor 2 Polyclonal Antibody

Cat No: HR1AP5811

For research use only

Overview

Product Name	Relaxin Receptor 2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	86kDa
GeneID?Human?	RXFP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RXFP2; GPR106; GREAT; LGR8; Relaxin receptor 2; G-protein coupled receptor 106; G-protein coupled receptor affecting testicular descent; Leucine-rich repeat-containing G-protein coupled receptor 8; Re
Background	relaxin/insulin like family peptide receptor 2(RXFP2) Homo sapiens This gene encodes a member of the GPCR (G protein-coupled, 7-transmembrane receptor) family. Mutations in this gene are associated with cryptorchidism. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Oct 2009],