

P2RY8 Polyclonal Antibody

Cat No: HR1AP4216

For research use only

Overview

Product Name	P2RY8 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	P2RY8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	P2RY8; P2Y purinoceptor 8; P2Y8
Background	<p>purinergic receptor P2Y8(P2RY8) Homo sapiens The protein encoded by this gene belongs to the family of G-protein coupled receptors, that are preferentially activated by adenosine and uridine nucleotides. This gene is moderately expressed in undifferentiated HL60 cells, and is located on both chromosomes X and Y. [provided by RefSeq, Jul 2008],</p>