

GPR115 Polyclonal Antibody

Cat No: HR1AP3560

For research use only

Overview

Product Name	GPR115 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	78kDa
GeneID?Human?	GPR115
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR115; PGR18; Probable G-protein coupled receptor 115; G-protein coupled receptor PGR18
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,