

GPR52 Polyclonal Antibody

Cat No: HR1AP3576

For research use only

Overview

Product Name	GPR52 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	GPR52
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR52; Probable G-protein coupled receptor 52
Background	G protein-coupled receptor 52(GPR52) Homo sapiens Members of the G protein-coupled receptor (GPR) family play important roles in signal transduction from the external environment to the inside of the cell.[supplied by OMIM, Jul 2002],