

GPR68 Polyclonal Antibody

Cat No: HR1AP8742

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | GPR68 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IF,ELISA |
| Species Reactivity | Human,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 41kDa |
| GeneID?Human? | GPR68 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | GPR68; OGR1; Ovarian cancer G-protein coupled receptor 1; OGR-1; G-protein coupled receptor 68; GPR12A; Sphingosylphosphorylcholine receptor |
| Background | similarity:Belongs to the G-protein coupled receptor 1 family., |