

GPR176 Polyclonal Antibody

Cat No: HR1AP3570

For research use only

Overview

Product Name	GPR176 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	GPR176
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR176; Probable G-protein coupled receptor 176; HB-954
Background	G protein-coupled receptor 176(GPR176) Homo sapiens Members of the G protein-coupled receptor family, such as GPR176, are cell surface receptors involved in responses to hormones, growth factors, and neurotransmitters (Hata et al., 1995 [PubMed 7893747]).[supplied by OMIM, Jul 2008],