

GPR137C Polyclonal Antibody

Cat No: HR1AP3564

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | GPR137C Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,IF,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 45kDa |
| GeneID?Human? | GPR137C |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | GPR137C; TM7SF1L2; Integral membrane protein GPR137C; Transmembrane 7 superfamily member 1-like 2 protein |
| Background | similarity:Belongs to the GPR137 family., |