

GPR15 Polyclonal Antibody

Cat No: HR1AP3565

For research use only

Overview

Product Name	GPR15 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	GPR15
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR15; G-protein coupled receptor 15; Brother of Bonzo; BoB
Background	G protein-coupled receptor 15(GPR15) Homo sapiens This gene encodes a G protein-coupled receptor that acts as a chemokine receptor for human immunodeficiency virus type 1 and 2. The encoded protein localizes to the cell membrane. [provided by RefSeq, Nov 2012],