

Axl (phospho Tyr691) Polyclonal Antibody

Cat No: HR1AP7881

For research use only

Overview

Product Name	Axl (phospho Tyr691) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	110 97kDa
GeneID?Human?	AXL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AXL; UFO; Tyrosine-protein kinase receptor UFO; AXL oncogene
Background	<p>AXL receptor tyrosine kinase (AXL) Homo sapiens The protein encoded by this gene is a member of the Tyro3-Axl-Mer (TAM) receptor tyrosine kinase subfamily. The encoded protein possesses an extracellular domain which is composed of two immunoglobulin-like motifs at the N-terminal, followed by two fibronectin type-III motifs. It transduces signals from the extracellular matrix into the cytoplasm by binding to the vitamin K-dependent protein growth arrest-specific 6 (Gas6). This gene may be involved in several cellular functions including growth, migration, aggregation and anti-inflammation in multiple cell types. Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, Jul 2013].</p>