

Glucagon Receptor Polyclonal Antibody

Cat No: HR1AP3544

For research use only

Overview

Product Name	Glucagon Receptor Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	GCGR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GCGR; Glucagon receptor; GL-R
Background	glucagon receptor(GCGR) Homo sapiens The protein encoded by this gene is a glucagon receptor that is important in controlling blood glucose levels. Defects in this gene are a cause of non-insulin-dependent diabetes mellitus (NIDDM).[provided by RefSeq, Jan 2010],