

## GPR101 Polyclonal Antibody

Cat No: HR1AP3555

For research use only

### Overview

Product Name	GPR101 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	GPR101
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR101; Probable G-protein coupled receptor 101
Background	G protein-coupled receptor 101(GPR101) Homo sapiens The protein encoded by this gene is an orphan G protein-coupled receptor of unknown function. The encoded protein is a member of a family of proteins that contain seven transmembrane domains and transduce extracellular signals through heterotrimeric G proteins. [provided by RefSeq, Sep 2011],