

MC5-R Polyclonal Antibody

Cat No: HR1AP3869

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | MC5-R Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IF,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 40kDa |
| GeneID?Human? | MC5R |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | MC5R; Melanocortin receptor 5; MC5-R; MC-2 |
| Background | melanocortin 5 receptor(MC5R) Homo sapiens This gene encodes a member of the seven-pass transmembrane G protein-coupled melanocortin receptor protein family that stimulate cAMP signal transduction. The encoded protein is a receptor for melanocyte-stimulating hormone and adrenocorticotrophic hormone and is suggested to play a role in sebum generation. [provided by RefSeq, Jun 2010], |