

FSHR Polyclonal Antibody

Cat No: HR1AP3489

For research use only

Overview

Product Name	FSHR Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	FSHR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FSHR; LGR1; Follicle-stimulating hormone receptor; FSH-R; Follitropin receptor
Background	follicle stimulating hormone receptor(FSHR) Homo sapiens The protein encoded by this gene belongs to family 1 of G-protein coupled receptors. It is the receptor for follicle stimulating hormone and functions in gonad development. Mutations in this gene cause ovarian dysgenesis type 1, and also ovarian hyperstimulation syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010],