

## SPINK6 Polyclonal Antibody

Cat No: HR1AP7250

For research use only

### Overview

Product Name	SPINK6 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	SPINK6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SPINK6; Serine protease inhibitor Kazal-type 6; Kallikrein inhibitor
Background	serine peptidase inhibitor, Kazal type 6(SPINK6) Homo sapiens The protein encoded by this gene is a Kazal-type serine protease inhibitor that acts on kallikrein-related peptidases in the skin. Two transcript variants the same protein have been found for this gene. [provided by RefSeq, Aug 2010],