

SH2D4A Polyclonal Antibody

Cat No: HR1AP8251

For research use only

Overview

Product Name	SH2D4A Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	SH2D4A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SH2D4A; SH2A; SH2 domain-containing protein 4A; Protein SH(2)A
Background	sequence caution:Translated as stop.,sequence caution:Translation N-terminally extended.,similarity:Contains 1 SH2 domain.,tissue specificity:Ubiquitously expressed. Aberrantly expressed in some cancers.,