

# Mic-1 Polyclonal Antibody

Cat No: HR1AP6810

For research use only

## Overview

Product Name	Mic-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	C18orf8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	C18orf8; MIC1; Uncharacterized protein C18orf8; Colon cancer-associated protein Mic1; Mic-1
Background	chromosome 18 open reading frame 8(C18orf8) Homo sapiens This gene encodes a colon cancer associated protein. [provided by RefSeq, Jan 2013],