

# Latherin Polyclonal Antibody

Cat No: HR1AP8168

For research use only

## Overview

Product Name	Latherin Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	LATH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LATH; BASE; Putative latherin; Breast cancer and salivary gland-expressed protein
Background	LATH, also named as BASE, belongs to the BPI/LBP/Plunc superfamily and Plunc family. It is the major protein in sweat, has surfactant properties.