

HLA-H Polyclonal Antibody

Cat No: HR1AP6853

For research use only

Overview

Product Name	HLA-H Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	HLA-H
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HLA-H; HLAH; Putative HLA class I histocompatibility antigen; alpha chain H; HLA-12.4; HLA-AR; MHC class I antigen H
Background	major histocompatibility complex, class I, H (pseudogene)(HLA-H) Homo sapiens This major histocompatibility complex gene represents a transcribed pseudogene, possibly derived from HLA-A. This gene displays extensive variation. [provided by RefSeq, May 2010],