

CD92 Polyclonal Antibody

Cat No: HR1AP6455

For research use only

Overview

Product Name	CD92 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	SLC44A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC44A1; CD92; CDW92; CTL1; Choline transporter-like protein 1; CDw92; Solute carrier family 44 member 1; CD antigen CD92
Background	function:Probable choline transporter. May be involved in membrane synthesis and myelin production.,similarity:Belongs to the CTL (choline transporter-like) family.,tissue specificity:Expressed in various cells of the hematopoietic system.,