

SH2D5 Polyclonal Antibody

Cat No: HR1AP4546

For research use only

Overview

Product Name	SH2D5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	SH2D5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SH2D5; SH2 domain-containing protein 5
Background	similarity:Contains 1 SH2 domain.,