

NEGR1 Polyclonal Antibody

Cat No: HR1AP4026

For research use only

Overview

Product Name	NEGR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	NEGR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NEGR1; IGLON4; Neuronal growth regulator 1; IgLON family member 4
Background	function:May be involved in cell-adhesion. May function as a trans-neural growth-promoting factor in regenerative axon sprouting in the mammalian brain.,similarity:Belongs to the immunoglobulin superfamily. IgLON family.,similarity:Contains 3 Ig-like C2-type (immunoglobulin-like) domains.,