

HLA-DO? Polyclonal Antibody

Cat No: HR1AP3648

For research use only

Overview

Product Name	HLA-DO? Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	HLA-DOB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HLA-DOB; HLA class II histocompatibility antigen; DO beta chain; MHC class II antigen DOB
Background	<p>major histocompatibility complex, class II, DO beta(HLA-DOB) Homo sapiens HLA-DOB belongs to the HLA class II beta chain paralogues. This class II molecule is a heterodimer consisting of an alpha (DOA) and a beta chain (DOB), both anchored in the membrane. It is located in intracellular vesicles. DO suppresses peptide loading of MHC class II molecules by inhibiting HLA-DM. Class II molecules are expressed in antigen presenting cells (APC: B lymphocytes, dendritic cells, macrophages). The beta chain is approximately 26-28 kDa and its gene contains 6 exons. Exon one encodes the leader peptide, exons 2 and 3 encode the two extracellular domains, exon 4 encodes the transmembrane domain and exon 5 encodes the cytoplasmic tail. [provided by RefSeq, Jul 2008],</p>