

Cleaved-C1s HC (R437) Polyclonal Antibody

Cat No: HR1AP2125

For research use only

Overview

Product Name	Cleaved-C1s HC (R437) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47+76kDa
GeneID?Human?	C1S
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	C1S; Complement C1s subcomponent; C1 esterase; Complement component 1 subcomponent s
Background	complement C1s(C1S) Homo sapiens This gene encodes a serine protease, which is a major constituent of the human complement subcomponent C1. C1s associates with two other complement components C1r and C1q in order to yield the first component of the serum complement system. Defects in this gene are the cause of selective C1s deficiency. [provided by RefSeq, Mar 2009],