

CD232 Polyclonal Antibody

Cat No: HR1AP5468

For research use only

Overview

Product Name	CD232 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	175kDa
GeneID?Human?	PLXNC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PLXNC1; VESPR; Plexin-C1; Virus-encoded semaphorin protein receptor; CD232
Background	<p>plexin C1(PLXNC1) Homo sapiens This gene encodes a member of the plexin family. Plexins are transmembrane receptors for semaphorins, a large family of proteins that regulate axon guidance, cell motility and migration, and the immune response. The encoded protein and its ligand regulate melanocyte adhesion, and viral semaphorins may modulate the immune response by binding to this receptor. The encoded protein may be a tumor suppressor protein for melanoma. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jan 2011],</p>