

# CD89 Polyclonal Antibody

Cat No: HR1AP5409

For research use only

## Overview

Product Name	CD89 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	FCAR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FCAR; CD89; Immunoglobulin alpha Fc receptor; IgA Fc receptor; CD89
Background	<p>Fc fragment of IgA receptor(FCAR) Homo sapiens This gene is a member of the immunoglobulin gene superfamily and encodes a receptor for the Fc region of IgA. The receptor is a transmembrane glycoprotein present on the surface of myeloid lineage cells such as neutrophils, monocytes, macrophages, and eosinophils, where it mediates immunologic responses to pathogens. It interacts with IgA-opsonized targets and triggers several immunologic defense processes, including phagocytosis, antibody-dependent cell-mediated cytotoxicity, and stimulation of the release of inflammatory mediators. Multiple alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008],</p>