

## C7 Polyclonal Antibody

Cat No: HR1AP9533

For research use only

### Overview

Product Name	C7 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	93kDa
GeneID?Human?	C7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	C7; Complement component C7
Background	complement C7(C7) Homo sapiens C7 is a component of the complement system. It participates in the formation of Membrane Attack Complex (MAC). People with C7 deficiency are prone to bacterial infection. [provided by RefSeq, Jul 2008],