

## NT5C1A Polyclonal Antibody

Cat No: HR1AP8877

For research use only

### Overview

Product Name	NT5C1A Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	NT5C1A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NT5C1A; Cytosolic 5'-nucleotidase 1A; cN1A; Cytosolic 5'-nucleotidase IA; cN-I; cN-IA
Background	5'-nucleotidase, cytosolic IA(NT5C1A) Homo sapiens Cytosolic nucleotidases, such as NT5C1A, dephosphorylate nucleoside monophosphates (Hunsucker et al., 2001 [PubMed 11133996]).[supplied by OMIM, Mar 2008],