

# NK-TR Polyclonal Antibody

Cat No: HR1AP7504

For research use only

## Overview

Product Name	NK-TR Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	NKTR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NKTR; NK-tumor recognition protein; NK-TR protein; Natural-killer cells cyclophilin-related protein
Background	natural killer cell triggering receptor(NKTR) Homo sapiens This gene encodes a membrane-anchored protein with a hydrophobic amino terminal domain and a cyclophilin-like PPIase domain. It is present on the surface of natural killer cells and facilitates their binding to targets. Its expression is regulated by IL2 activation of the cells. [provided by RefSeq, Jul 2008],