

TCR ? Polyclonal Antibody

Cat No: HR1AP4687

For research use only

Overview

Product Name	TCR ? Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	TRBC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TRBC1; T-cell receptor beta-1 chain C region
Background	similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,