

CD72 Polyclonal Antibody

Cat No: HR1AP9499

For research use only

Overview

Product Name	CD72 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	CD72
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD72; B-cell differentiation antigen CD72; Lyb-2; CD72
Background	function:Plays a role in B-cell proliferation and differentiation. Associates with CD5.,online information:CD72,similarity:Contains 1 C-type lectin domain.,subunit:Homodimer; disulfide-linked.,tissue specificity:Pre-B-cells and B-cells but not terminally differentiated plasma cells.,