

MBL-C Polyclonal Antibody

Cat No: HR1AP9649

For research use only

Overview

Product Name	MBL-C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	MBL2 COLEC1 MBL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MBL2 COLEC1 MBL
Background	<p>mannose binding lectin 2(MBL2) Homo sapiens This gene encodes the soluble mannose-binding lectin or mannose-binding protein found in serum. The protein encoded belongs to the collectin family and is an important element in the innate immune system. The protein recognizes mannose and N-acetylglucosamine on many microorganisms, and is capable of activating the classical complement pathway. Deficiencies of this gene have been associated with susceptibility to autoimmune and infectious diseases. [provided by RefSeq, Jul 2008].</p>