

## NET-5 Polyclonal Antibody

Cat No: HR1AP4031

For research use only

### Overview

Product Name	NET-5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	TSPAN9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TSPAN9; NET5; Tetraspanin-9; Tspan-9; Tetraspan NET-5
Background	<p>tetraspanin 9(TSPAN9) Homo sapiens The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. Alternatively spliced transcripts encoding the same protein have been identified. [provided by RefSeq, Nov 2009],</p>