

CD5 (phospho Tyr453) Polyclonal Antibody

Cat No: HR1AP9100

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | CD5 (phospho Tyr453) Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,IF,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 40kDa |
| GeneID?Human? | CD5 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | CD5; LEU1; T-cell surface glycoprotein CD5; Lymphocyte antigen T1/Leu-1; CD antigen CD5 |
| Background | function:May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.,similarity:Contains 3 SRCR domains., |