

# A1BG Polyclonal Antibody

Cat No: HR1AP4898

For research use only

## Overview

Product Name	A1BG Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	A1BG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	A1BG; Alpha-1B-glycoprotein; Alpha-1-B glycoprotein
Background	alpha-1-B glycoprotein(A1BG) Homo sapiens The protein encoded by this gene is a plasma glycoprotein of unknown function. The protein shows sequence similarity to the variable regions of some immunoglobulin supergene family member proteins. [provided by RefSeq, Jul 2008],