

Ox40L Polyclonal Antibody

Cat No: HR1AP5461

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | Ox40L Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 21kDa |
| GeneID?Human? | TNFSF4 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | TNFSF4; TXGP1; Tumor necrosis factor ligand superfamily member 4; Glycoprotein Gp34; OX40 ligand; OX40L; TAX transcriptionally-activated glycoprotein 1; CD252 |
| Background | tumor necrosis factor superfamily member 4(TNFSF4) Homo sapiens This gene encodes a cytokine of the tumor necrosis factor (TNF) ligand family. The encoded protein functions in T cell antigen-presenting cell (APC) interactions and mediates adhesion of activated T cells to endothelial cells. Polymorphisms in this gene have been associated with Sjogren's syndrome and systemic lupus erythematosus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014], |