

CD34 Polyclonal Antibody

Cat No: HR1AP3018

For research use only

Overview

Product Name	CD34 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	CD34
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD34; Hematopoietic progenitor cell antigen CD34; CD antigen CD34
Background	CD34 molecule(CD34) Homo sapiens The protein encoded by this gene may play a role in the attachment of stem cells to the bone marrow extracellular matrix or to stromal cells. This single-pass membrane protein is highly glycosylated and phosphorylated by protein kinase C. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011],