

# LMTK3 Polyclonal Antibody

Cat No: HR1AP5759

For research use only

## Overview

|                       |   |
|-----------------------|---|
| Product Name          | LMTK3 Polyclonal Antibody   |
| Source                | Rabbit  |
| Applications          | WB,IHC-p,ELISA  |
| Species Reactivity    | Human,Mouse   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 40kDa   |
| GeneID?Human?         | LMTK3   |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | LMTK3; KIAA1883; TYKLM3; Serine/threonine-protein kinase LMTK3; Lemur tyrosine kinase 3   |
| Background            | catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Tyr protein kinase family.,similarity:Contains 1 protein kinase domain., |