

ADCK4 Polyclonal Antibody

Cat No: HR1AP2713

For research use only

Overview

Product Name	ADCK4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	ADCK4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ADCK4; Uncharacterized aarF domain-containing protein kinase 4
Background	<p>coenzyme Q8B(COQ8B) Homo sapiens This gene encodes a protein with two copies of a domain found in protein kinases. The encoded protein has a complete protein kinase catalytic domain, and a truncated domain that contains only the active and binding sites of the protein kinase domain, however, it is not known whether the protein has any kinase activity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011],</p>