

RIP140 (Acetyl Lys158) Polyclonal Antibody

Cat No: HR1AP9545

For research use only

Overview

Product Name	RIP140 (Acetyl Lys158) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	NRIP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NRIP1; Nuclear receptor-interacting protein 1; Nuclear factor RIP140; Receptor-interacting protein 140
Background	nuclear receptor interacting protein 1(NRIP1) Homo sapiens Nuclear receptor interacting protein 1 (NRIP1) is a nuclear protein that specifically interacts with the hormone-dependent activation domain AF2 of nuclear receptors. Also known as RIP140, this protein modulates transcriptional activity of the estrogen receptor. [provided by RefSeq, Jul 2008],