

Lyn (phospho Tyr508) Polyclonal Antibody

Cat No: HR1AP2585

For research use only

Overview

Product Name	Lyn (phospho Tyr508) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53,56kDa
GeneID?Human?	LYN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LYN; JTK8; Tyrosine-protein kinase Lyn; Lck/Yes-related novel protein tyrosine kinase; V-yes-1 Yamaguchi sarcoma viral related oncogene homolog; p53Lyn; p56Lyn
Background	LYN proto-oncogene, Src family tyrosine kinase(LYN) Homo sapiens This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011],