

QSK Polyclonal Antibody

Cat No: HR1AP6447

For research use only

Overview

Product Name	QSK Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	140kDa
GeneID?Human?	SIK3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SIK3; KIAA0999; QSK; L19; Serine/threonine-protein kinase SIK3; Salt-inducible kinase 3; SIK-3; Serine/threonine-protein kinase QSK
Background	similarity:Belongs to the protein kinase superfamily.,