

## DAPK3 (phospho Thr265) Polyclonal Antibody

Cat No: HR1AP6078

For research use only

### Overview

Product Name	DAPK3 (phospho Thr265) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	DAPK3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DAPK3; ZIPK; Death-associated protein kinase 3; DAP kinase 3; DAP-like kinase; Dlk; MYPT1 kinase; Zipper-interacting protein kinase; ZIP-kinase
Background	death associated protein kinase 3(DAPK3) Homo sapiens Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis. [provided by RefSeq, Jul 2008],