

RIOK1 Polyclonal Antibody

Cat No: HR1AP8816

For research use only

Overview

Product Name	RIOK1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RIOK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RIOK1; Serine/threonine-protein kinase RIO1; RIO kinase 1
Background	RIO kinase 1(RIOK1) Homo sapiens This gene includes two alternatively spliced transcript variants, which encode different isoforms. The function of this gene has not been determined. [provided by RefSeq, Jul 2008],