

RPS6KL1 Polyclonal Antibody

Cat No: HR1AP8812

For research use only

Overview

Product Name	RPS6KL1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RPS6KL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPS6KL1; Ribosomal protein S6 kinase-like 1
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family. S6 kinase subfamily.,similarity:Contains 1 MIT domain.,similarity:Contains 1 protein kinase domain.,