

PRKY Polyclonal Antibody

Cat No: HR1AP7953

For research use only

Overview

Product Name	PRKY Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	PRKY
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PRKY; Serine/threonine-protein kinase PRKY
Background	<p>protein kinase, Y-linked, pseudogene(PRKY) Homo sapiens This gene is similar to the protein kinase, X-linked gene in the pseudoautosomal region of the X chromosome. The gene is classified as a transcribed pseudogene because it has lost a coding exon that results in all transcripts being candidates for nonsense-mediated decay (NMD) and unlikely to express a protein. Abnormal recombination between this gene and a related gene on chromosome X is a frequent cause of XX males and XY females. [provided by RefSeq, Jul 2010].</p>