

## Pim-2 Polyclonal Antibody

Cat No: HR1AP5697

For research use only

### Overview

Product Name	Pim-2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PIM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PIM2; Serine/threonine-protein kinase pim-2; Pim-2h
Background	Pim-2 proto-oncogene, serine/threonine kinase(PIM2) Homo sapiens This gene encodes a protooncogene that acts as a serine/threonine protein kinase. Studies determined the encoded protein functions to prevent apoptosis and to promote cell survival.[provided by RefSeq, Nov 2009],