

PHACTR4 Polyclonal Antibody

Cat No: HR1AP8350

For research use only

Overview

Product Name	PHACTR4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	PHACTR4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PHACTR4; Phosphatase and actin regulator 4
Background	<p>phosphatase and actin regulator 4(PHACTR4) Homo sapiens This gene encodes a member of the phosphatase and actin regulator (PHACTR) family. Other PHACTR family members have been shown to inhibit protein phosphatase 1 (PP1) activity, and the homolog of this gene in the mouse has been shown to interact with actin and PP1. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],</p>