

MAST4 Polyclonal Antibody

Cat No: HR1AP7112

For research use only

Overview

Product Name	MAST4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	MAST4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MAST4; KIAA0303; Microtubule-associated serine/threonine-protein kinase 4
Background	microtubule associated serine/threonine kinase family member 4(MAST4) Homo sapiens This gene encodes a member of the microtubule-associated serine/threonine protein kinases. The proteins in this family contain a domain that gives the kinase the ability to determine its own scaffold to control the effects of their kinase activities. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2014],