

Jamip2 Polyclonal Antibody

Cat No: HR1AP9215

For research use only

Overview

Product Name	Jamip2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	JAKMIP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	JAKMIP2; JAMIP2; KIAA0555; NECC1; Janus kinase and microtubule-interacting protein 2; CTCL tumor antigen HD-CL-04; Neuroendocrine long coiled-coil protein 1
Background	janus kinase and microtubule interacting protein 2(JAKMIP2) Homo sapiens The protein encoded by this gene is reported to be a component of the Golgi matrix. It may act as a golgin protein by negatively regulating transit of secretory cargo and by acting as a structural scaffold of the Golgi. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012],